

Laboratory True Calibration Services, 68, Silicon Industrial Hub, Moraiya, Ahmedabad, Gujarat

Accreditation Standard ISO/IEC 17025: 2005

Certificate Number CC-2494 Page 1 of 5

Validity 16.08.2018 to 15.08.2020 Last Amended on -

Sl.	Quantity Measured / Instrument	Range/Frequency	*Calibration Measurement Capability (\pm)	Remarks
<u>MECHANICAL CALIBRATION</u>				
I.	WEIGHTS			
1.	Standard Weights ^s	1 mg 2 mg 5 mg 10 mg 20 mg 50 mg 100 mg 200 mg 500 mg 1 g 2 g 5 g 10 g 20 g 50 g	0.02 mg 0.02 mg 0.02 mg 0.02 mg 0.03 mg 0.04 mg 0.05 mg 0.05 mg 0.05 mg 0.05 mg 0.10 mg 0.10 mg 0.10 mg 0.10 mg 0.10 mg 0.10 mg	Using E2 Standard Weights & Semi-Micro Balance (Readability: 0.01 mg) Calibration of Weights of Class F2 Accuracy and Coarser as per OIML R-111
		100 g 200 g	0.4 mg 0.4 mg	Using E2 Standard Weights & Analytical Balance (Readability: 0.1 mg) Calibration of Weights of Class F2 Accuracy and Coarser as per OIML R-111
II.	VOLUME			
1.	Micro Pipette ^s	10 μ l to 100 μ l @ 27 °C >100 μ l to 1000 μ l @ 27 °C >1000 μ l to 5000 μ l @ 27 °C	0.26 μ l 0.50 μ l 7.10 μ l	Using Weighing Balance of 80 g Capacity and 0.01 mg Readability by Gravimetric Method Based on ISO 8655 Part 6

Ram Ashray
Convenor

Avijit Das
Program Manager

Laboratory True Calibration Services, 68, Silicon Industrial Hub, Moraiya, Ahmedabad, Gujarat

Accreditation Standard ISO/IEC 17025: 2005

Certificate Number CC-2494 **Page** 2 of 5

Validity 16.08.2018 to 15.08.2020 **Last Amended on** -

Sl.	Quantity Measured / Instrument	Range/Frequency	*Calibration Measurement Capability (\pm)	Remarks
		>5000 μ l to 10000 μ l @ 27 °C	7.10 μ l	
2.	Glassware Like Pipettes, Burettes, Measuring Cylinder, Volumetric Flask etc ^{\$}	1 ml to 10 ml @27 °C >10 ml to 100 ml @27 °C	0.12 ml 0.22 ml	Using Weighing Balance of 220 g Capacity and 0.1 mg Readability by Gravimetric Method Based on IS/ISO 4787
III.	WEIGHING SCALE & BALANCE			
1.	Electronic Weighing Balances with readability d=0.01 mg [*]	Maximum capacity Up to 80 g	0.05 mg	Using E2 class weights Calibration of electronic weighing balance and Comparator of Class I and coarser As per OIML R-76-1
2.	Electronic Weighing Balances with readability d=0.1 mg [*]	Maximum capacity Up to 200 g	0.11 mg	Using E2 class weights Calibration of electronic weighing balance and Comparator of Class I and coarser As per OIML R-76-1
3.	Electronic Weighing Balances with Readability d=0.1 g [*]	Maximum Capacity Up to 5000 g	0.1 g	Using F1 class weights Calibration of electronic weighing balance and Comparator of Class III and coarser As per OIML R-76-1
IV.	ACCELERATION AND SPEED			
1.	Centrifuge, Shaker, Rotor [*]	100 RPM to 15000 RPM	6.8 RPM to 59.0 RPM	Using Digital Tachometer by Comparison Method

Ram Ashray
Convenor

Avijit Das
Program Manager

Laboratory True Calibration Services, 68, Silicon Industrial Hub, Moraiya,
 Ahmedabad, Gujarat
Accreditation Standard ISO/IEC 17025: 2005
Certificate Number CC-2494 **Page** 3 of 5
Validity 16.08.2018 to 15.08.2020 **Last Amended on** -

Sl.	Quantity Measured / Instrument	Range/Frequency	*Calibration Measurement Capability (\pm)	Remarks
V.	PRESSURE INDICATING DEVICES			
1.	Digital & Analog Pressure Gauge/ Transmitter/Indicator ^s	0 to 20 bar 0 to 700 bar	1.05 bar 3.12 bar	Using Digital Pressure Calibrator by Comparison method as per DKD-R 6-1
2.	Digital & Analog Pressure Gauge/ Transmitter/Indicator ^s	0 to 700 bar	2.68 bar	Using Digital Pressure Calibrator by Comparison method as per DKD-R 6-1

Ram Ashray
Convenor

Avijit Das
Program Manager

Laboratory True Calibration Services, 68, Silicon Industrial Hub, Moraiya, Ahmedabad, Gujarat

Accreditation Standard ISO/IEC 17025: 2005

Certificate Number CC-2494 **Page** 4 of 5

Validity 16.08.2018 to 15.08.2020 **Last Amended on** -

Sl.	Quantity Measured / Instrument	Range/Frequency	*Calibration Measurement Capability (\pm)	Remarks
<u>THERMAL CALIBRATION</u>				
1.	TEMPERATURE			
1.	Liquid in Glass Thermometer [§]	50 °C to 150 °C	0.39 °C	Using Precision Temperature Scanner & SSPRT(Fluke 1586 & 5609) & Oil bath by Comparison Method
2.	RTD, Thermocouple, Thermistor, Digital Thermometer, Temperature Controller/ Indicator with Sensor, Temperature Transmitter, Temperature Gauge [§]	(-) 80 °C to (-) 20 °C (-) 20 °C to 50 °C 50 °C to 500 °C	0.47 °C 0.47 °C 0.53 °C	Using Precision Temperature Scanner & SSPRT(Fluke 1586 & 5609) with Liquid Bath, Oil Bath and Dry Block by Comparison Method
3.	RTD, Thermocouple, Thermistor, Digital Thermometer, Temperature Controller/Indicator with Sensor, Temperature Transmitter, Temperature Gauge [*]	(-) 25 °C to 50 °C 50 °C to 500 °C	0.47 °C 0.59 °C	Using Precision Thermometer, Multifunction Calibrator, Thermocouple with Dry Block by Comparison Method
4.	Deep Freezer, Freezer Incubator (for Non-Medical Applications), Oven, Autoclave (Single Position) [*]	(-) 80 °C to 10 °C 10 °C to 100 °C 100 °C to 300 °C	1.20 °C 0.40 °C 0.70 °C	Using Precision Thermometer, Multi-Function Calibrator (Beamex MC2-TE), Process Data Logger, RTD Sensor by Single

Ram Ashray
Convenor

Avijit Das
Program Manager

Laboratory True Calibration Services, 68, Silicon Industrial Hub, Moraiya, Ahmedabad, Gujarat
Accreditation Standard ISO/IEC 17025: 2005
Certificate Number CC-2494 **Page** 5 of 5
Validity 16.08.2018 to 15.08.2020 **Last Amended on** -

Sl.	Quantity Measured / Instrument	Range/Frequency	*Calibration Measurement Capability (\pm)	Remarks
				Position Calibration Method
5.	Oven, Muffle Furnace*	500 °C to 1200 °C	3.0 °C	Using Multi-function Calibrator, S type Thermocouple by Single Position Calibration Method
6.	Deep Freezer, Freezer Incubator (for Non-Medical Applications), Oven, Autoclave (for Non-Medical Applications) Cold Room, Clean Room, Environmental Chamber (Multi Position)*	(-) 80 °C to 250 °C	2.73 °C	Using 16 Channel Data Logger (Masibus 85XX) by Multi Position Method
II.	SPECIFIC HEAT AND HUMIDITY			
1.	Thermo Hygrometer, Humidity Data Logger/ Transmitter/Indicator [§]	30 % RH to 95 % RH 10 °C to 50 °C	2.90 % RH 0.57 °C	Using Temperature & Humidity Indicator (Rotronic/Polltech) with Humidity Chamber by Comparison Method
2.	Humidity Chamber, Stability Chamber, Environmental Chamber*	30 % RH to 95 % RH	4.74 % RH @ 25 °C	Using Temperature & Humidity Datalogger 10 Nos. by Multi Position Method

* Measurement Capability is expressed as an uncertainty (\pm) at a confidence probability of 95%

[§]Only in Permanent Laboratory

^{*}Only for Site Calibration

Ram Ashray
 Convenor

Avijit Das
 Program Manager